

# ADVANCE BEAUTY

ADVANCE **E**BEAUTY



Rovaltain, parc 45 parallèle  
5 Rue Olivier de Serres  
26958 Valence Cedex 9  
FRANCE

Tél : +33 (0)4 75 48 17 54  
Fax : +33 (0)4 75 47 74 93

Info@advance-beauty.com  
www.advance-beauty.com



Création : MACOMAMOL.FR

MADE IN FRANCE

## CAVITA-LYSE **Sequential**

Cavitation is now the reference for efficient slimming localized treatments and provides a high profitability practice.

Thanks to its patented Sequential technology, the action processes are maximized and the results are sure to come from the first session.



ADVANCE **E**BEAUTY



# WHAT IS CAVITATION?

## PHYSICAL PRINCIPLE OF CAVITATION

Ultrasound waves are elastic waves that come from sequences of compression and depression of a piezoelectric material, emitted through a transducer in contact with the skin. When a liquid is subjected to an ultrasonic mechanical vibration having right parameters, it appears a set of pressions and depressions creating «micro-cavities», that means small gaz bubbles. These cavities grow and implode leading locally very high temperature and internal pressure of thousandfolds the atmospheric pressure.

This phenomenon is called CAVITATION.

## BIOLOGICAL EFFECTS OF CAVITATION

The implosion of small bubbles into the adipose tissue release, at the molecular scale, a very powerful energy that creates impact waves in the interstitial liquids. This energy is able to break surrounding solid particles creating the lysis of adipose tissue, a destruction of cells localized in fat tissue.

### CAVITATION IS A KIND OF NON INVASIVE LIPOPLASTY.

CAVITATION probe is used to get powerful energy on the localized fat. Its titanium resonator is fully biocompatible and is a perfect vibrations transducer.

CAVITA-LYSE  Sequential

CAVITA-LYSE Sequential get remarkable results for reshaping and slimming programs: reduction of adipose tissue, cellulite, orange peel skin, firming...

Dedicated to beauty salons and slimming centers, and coming from a new technique totally controlled, CAVITA-LYSE Sequential is very easy to program and use.

This NON INVASIVE fat reduction process needs only sessions of less than 30 minutes.

Sequential cavitation provides a localized but not focalized action (due to its low frequency only), so we apply a safe treatment and get the most efficient results.

CAVITA-LYSE Sequential is the first cavitation device all over made in France.

Thanks to the latest technology of its titanium probe, CAVITA-LYSE Sequential insure a safe emission.

## WHY THE SEQUENTIAL TECHNOLOGY?

Stimulated by mechanical waves the physiological processes which act on the skin will be “reactivated”, improving micro-circulations and implementing fat destocking.

The Sequential mode principle, patented by Advance Beauty, is based on a chain of ultrasonic sequences, different and non-repetitive. So this “reactivation” remains constant throughout the treatment to increase the results and guarantee maximum efficiency for the full duration of the session, without habituation.

## PROVEN RESULTS

AFTER 3 SESSIONS ONLY.

A significant loss in circumference:

- 3.88 

ON THE ABDOMEN

- 4.35 

ON THE HIPS

\* Clinical study on 28 women



## TECHNICAL CHARACTERISTICS

Electricity supply	From 100V to 240 V AC / 50-60 Hz
Dimensions (cm)	30 (h) x 49(l) x 28(p) cm
3 modes	Sequential (exclusivity Advance Beauty), Firming & Shapping
Weight	6 Kg

## PROBE SPECIFICATIONS

CAVITATION Probe Frequency	40 KHz
Maximum real acoustic power	25 W
Active treatment area	20 cm²
Weight	700 g

