

DOCTOR LIFE[®]

4MAX

Strong but Silent
Homecare Presso-system


Improve your performance



What's advantage?

4MAX
Strong but Silent
Homecare Presso-system





Features

4MAX
Strong but Silent
Homecare Presso-system



DOCTOR LIFE®

4 modes
setting
(Massage/Circulation/Reverse/Full)

Features

4MAX
Strong but Silent
Homecare Presso-system



4MAX QUICK GUIDE

! SAFE USE TIP

Check your body condition at the lowest pressure for 10 minutes if you are a beginner, and then increase the setting gradually.


PRESSURE



Low : 1~2 Level
Mid : 3~6 Level
High : 7~10 Level

• Find a comfortable pressure according to your body condition.

SKIP



• You can skip a chamber(s) you don't want.

MODE



WAVE




CIRCULATION



UP/DOWN



FULL



• [UP/DOWN mode] operates the [Reverse WAVE mode] after the [CIRCULATION mode] twice.

Order package?

4MAX
Strong but Silent
Homecare Presso-system

Leg set (components)



DEVICE



LEG SLEEVE
(L, XL, XXL size)



DOUBLE TUBING




USER MANUAL

Leg sleeve advantage

4MAX
Strong but Silent
Homecare Presso-system

Comfortable
sole(foot)
during
compression.



One-body
hose cover

Foot
chamber

Casual pattern



Optional accessory?

4MAX
Strong but Silent
Homecare Presso-system



ARM SLEEVE



LEG SLEEVE



WAIST SLEEVE




RECOVERY BAG



STRETCHING MAT 4



CENTER BODY SLEEVE



Thank you!